



## Sabina™ 200 Lift



### For simple, reliable raising and transfer

The Sabina 200 Lift is used to assist people who have difficulty rising from a sitting to a standing position, for example from a wheelchair, chair, bed or when using the toilet.

The lift mast has three settings to match the height of the individual patient. The carefully designed accessories ensure that rising and stand related settings can be adapted to suit the needs of the individual patient. This improves the conditions for the safest and most comfortable raising possible.

## The best way of raising and transferring your patient

The Sabina™ lift is the result of comprehensive study and careful analysis of ergonomic principles and body movement patterns. The Sabina lift's lifting function and anchor points close to the patient's center of gravity ensure the most effective way of raising and transfer.

### Accessories for adaptation to the individual

There are three types of sit-to-stand vests available; select the one which suits your patient best. If the patient requires a little extra help to start the lift, you may wish to use a Seatstrap. The Sabina SideSupport is another alternative for patients with a tendency to lean to the side.

### Accessories

2027101	Sabina SideSupport
2027011	Heel support Sabina
2027006	Seatstrap Slingbar
3591115	Sabina Seatstrap



The emergency stop can be easily reached, but its position on the side reduces the risk of it being pressed accidentally. The electrical emergency lowering control is also easily reached.

### ■ TECHNICAL DATA

Max. Load:	350 lbs (160 kg)
Lifting Speed:	0.9 in./sec (20 mm/sec) without load
Batteries:	2 pcs. 12 V 2.9 Ah valve-regulated closed lead acc. so called gel batteries
Battery Charger:	External charger, 100-240 V AC 50-60 Hz, max 220 mA
Front Wheels:	Standard 3" (75 mm) twin wheel
Rear Wheels:	Standard 3" (75 mm) twin wheel with brakes
Emergency Lowering:	Electrical
Protection Class:	IPX4
Material:	Powder-coated Steel
Weight:	92 lbs (42 kg)
Product Number:	2020020 Sabina 200



Raising and base-width adjustment are easily controlled by the caregiver via the hand control.



The Sabina 200 lift has been tested by an accredited testing institute and complies with the requirements for Medical Device Directive Class I products (MDD 93/42/EEC).

The Sabina 200 lift complies with the requirements for IEC 60601-1, IEC 60601-1-2, EN ISO 10535, UL-60601-1 and CAN/CSA C22.2 No.601.1.

Design and Quality by Liko in Sweden.

Liko is quality-certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485.

Liko is also certified in accordance with environmental standard ISO 14001.

